Exploring with Tech Trek Microsco Student Guide	opes Name				
Look at two or three specimens using the different magnifying devices available. Then choose one specimen. Record your observations. (Size, color, shape, apparent texture, other descriptions of what you see. Stick to the facts. Comparative information may help you "recall" your observations, but are not the same as an actual observation. Describe "what you see" using each device. Record your data by drawing descriptive pictures in the box or circle provided.					
PSEM (Personal Scanning Electron Microscope) or SEM					
	Observations				
Specimen #1					
_					
_					
Magnification					
TYPES OF LIGHT MICROSCOPES					
Intel Play (BLUE) computer aided microscope Observations					
	Specimen #2				

Magnification_____

Compound Binocular Light Microscope or Medical Type					
Label information from glass slide					
Observations	Magnification (x)				
Specimen #3 Label information					
Dissection Microscope (65X for all) CSI TYPE Observations					
	Specimen #4				

Exploring Tech Trek	Microscopes	Name	
Student Guide			

Use the materials available to create a presentation, which includes a chart, to share your findings with the class.

Your presentation should include, but is not limited to:

- Photos or pictures of specimens
- Advantages/disadvantages of each magnifying device
- Similarities and differences in specimens when viewed under the different magnifiers.

Your presentation may also include:

- Information about TEM, SEM, ESEM, and AFM.
- Uses for each different type of microscope.
- Jobs associated with the different types of microscopes.

Your presentation should be:

- Mounted on poster-board
- Clearly labeled
- Grammatically accurate
- Visually pleasing
- An accurate display of information

Your presentation may be:

• Presented with multimedia (PowerPoint, etc.)

Multimedia presentations should be:

- Clearly labeled
- Grammatically accurate
- Visually pleasing
- An accurate display of information